








Slings

Human Care Lifts has a wide range of slings for a variety of users and their ever-changing needs. The slings are constantly being developed in conjunction with customers and experts in the field. Documented testing is the first and most important step for attaining the optimal lift. We are happy to help with advising, demonstrations, training and testing. You can find our representative for your area on our website, where you can also find up-to-date information, print out instructions, product information and more. www.humancarelift.com

Colour-codes:	Junior Small (JS)		Junior Medium (JM)		Junior Large (JL)		
	Small (S)		Medium (M)		Large (L)		Extra Large (XL) 

Transfer slings



Full Body Sling Plus 25010



Sit Sling 25030

Full Body Sling Plus 25010 (S, M, L, XL)

Approved for two- and three-point suspension

The sling provides a reclined sitting position. Together with the reinforced, tailored neck support, it gives optimum support throughout the transfer process. Works well for people with reduced neck mobility and poor body stability.

Sit Sling 25030 (JS, JM, S, M, L, XL)

Approved for two- and three-point suspension

Option: Human Care Breast Belt 25145/25146

A sling which suits most people. It has a reinforced, half-height back support and gives an upright sitting position with excellent support for the entire back. Its simplicity of use makes this one of our most popular slings.

Toileting slings



Toileting Sling Plus 25080



Toileting Sling Plus net 25090

Toileting Sling Plus 25080 (JS, JM, S, M, L, XL)

Toileting Sling Plus net 25085 (S, M, L, XL)

Approved for two- and three-point suspension

Option: Human Care Breast Belt 25145/25146

Simple and easy-to-use slings which allow for dressing and undressing associated with lifts to and from the toilet. The reinforced back support and strong padding under the arms provide comfortable support. The Breast Belt option gives enhanced safety and comfort.

Toileting Sling High net 25090 (S, M, L)

Approved for two- and three-point suspension

Option: Human Care Breast Belt 25145/25146

The sling has the same cut in the lower area as the Toileting Sling Plus. The neck support is height adjustable, allowing it to work perfectly regardless of the height of the user. The Breast Belt option gives enhanced safety and comfort.

Amputee/ safety



Amputee Sling 25040



Bath Sling 25035

Bath

Amputee Sling 25040 (S, M, L, XL)

Approved for two- and three-point suspension

Option: Human Care Breast Belt 25145/25146

Option: Human Care Lift spreader bar 21126/50268 for four-point suspension

A sling suitable for users with many different needs. For amputees the leg sections should be overlapped under both legs and inserted through the special loops on the inside and outside of the sling (see picture). The cut of the sling also makes it suitable for users who are at greater risk of sliding out, as well as with users who have hip problems, which calls for an open hip angle. The sling's leg sections can then be overlapped or pulled between the legs.

Bath Sling 25035 (S, M, L, XL)

Approved for two- and three-point suspension

A sling specially suited to bathing and showering. The material is a polyester net. The back is reinforced with two plastic slats to give good, stable support. If neck support is needed, we recommend Net Sling High 25015.

Net slings



Net Sling High 25015



Soft Sling 25050

Net Sling High 25015
Soft Sling High 25020
(JS, JM, S, M, L, XL)

Net Sling 25045
Soft Sling 25050

Approved for two- and three-point suspension

Soft and supple slings without handles and stiffeners which give the possibility of leaving them on when sitting in a chair. The high model gives a reclined sitting position with comfortable support for the neck, while the low model stops at the shoulders and gives a more upright sitting position.

25015/25045 are made of a low-friction net suitable for narrow wheelchairs and specially-shaped hammocks, as well as for bathing and showering.

25020/25050 are made of a high-friction, cotton-like polyester to minimize the risk of the user sliding out.

Comfort slings



Comfort Sling 25060



Comfort Sling High 25065

Comfort Sling 25060 (JM, S, M, L, XL)
Comfort Sling High 25065 (S, M, L, XL)

Approved for two- and four-point suspension

Option: Human Care Lift spreader bar 21126/50268

Soft and supple slings without handles and stiffeners. The comfort slings have no separate leg support and are put on with the user lying down, and then remain on the user when sitting. The slings give a equal weight distribution over a large area. Comfort Sling High has a shaped neck support.

Active



Activity Sling 25100



Multisling 25105

Activity Sling 25100 (JS, JM, JL)

Approved for two-point suspension

A junior sling which makes it possible for the user to get up and move about according to his or her own abilities. It has a cup-shaped neck support and elastic bands on the sides for optimal flexibility.

Multisling 25105 (S, M, L)

Approved for two-, three- and four-point suspension

The sling consists of five different parts. Use the breast belt as support when doing walking and balancing exercises, going to the toilet or dressing/undressing. Combine the Breast Belt with the Leg Support for additional support and security. When exercising with patients with considerable functional disorders, eg stroke-induced, the sling can be supplemented with the Hip Belt (25107). A breast band and shoulder padding are also included with the Multisling.

Tailor-made slings

All models have unique attributes to allow them to work in different lifting situations. If you need any special changes or adjustments, please contact our representative in your area.

HUMAN
CARE
Lift